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IMPORTANT : INFORMATION NOTE

Dear customer,

First of all we would like to thank you for choosing a Stago analyzer, and most of all we would like to remind you the importance that Diagnostica Stago gives to the quality of its products.

Each product undergoing an extensive validation before commercialisation.

As you know, our analyzers function with their associated consumables, particularly with the cuvettes in which the reaction between patient's plasma and the reagent takes place.

Quality and reliability of the result depend on different parameters like the technical specifications of these cuvettes. That's why we apply a very strict quality validation of the whole process of their production.

This validation process includes among other steps:

- a control of the proper functioning of the cuvette in the "Stago system" which includes the analyzer, the reagents and the consumables.

The control of our production process and quality control process allow to ensure the compatibility of all these elements and to provide the customers a guarantee of the quality and conformity to our specifications, on all our products.

We have been informed that some customers are using on STart® 4 analyzers cuvettes strips that are not manufactured, not validated and thus not guaranteed by Diagnostica Stago.

We would like to draw your attention on the risks caused by this usage for your patients and consequently for your laboratory.

The compatibility of these cuvettes strips with the system is not validated, thus their usage cancels the application of the Stago guarantee on your Stago system. So Stago declines any responsibility in case of dysfunction or erroneous diagnosis caused by the use of these cuvettes.

Moreover Stago has realized comparison testing between Stago cuvettes and non validated cuvettes coming from the Moroccan market. **These tests have shown that non validated cuvettes:**

- **do not fulfil Stago acceptance criteria**
- **give results that may lead to error in diagnostic (see data on page 2)**

In the case your laboratory is offered cuvettes not validated by Diagnostica Stago, we ask you to inform your official distributor so that we can take actions to ensure the quality of the products you are using to conduct your haemostasis analysis.

We remain at your disposal for any additional information you may need.

Best regards,

Nayla Ayoub
Product Manager
International Distribution Network

A handwritten signature in blue ink, appearing to read "Nayla Ayoub".

Comparison:

Tests results obtained with STAGO cuvettes and duplicate (non validated) cuvettes

Repeatability tests:

n = 24

Parameter: Prothrombin Time (PT)

Level: Normal

STAGO acceptance criteria: **CV < 1.5%**

	CV obtained with STAGO cuvettes	CV obtained with non validated cuvettes
STart 4 analyzer n°1	1,21	2,21
STart 4 analyzer n°2	0,83	1,64


n = 24

Parameter: Fibrinogen

Level : Pathhologic

STAGO acceptance criteria: **CV < 5%**

	CV obtained with STAGO cuvettes	CV obtained with non validated cuvettes
STart 4 analyzer n°1	2,99	16,14
STart 4 analyzer n°2	2,71	21,51

As you can observe from the graph, the same sample can vary from 0,82 to 1,88 g/L or 1,07 to 1,82 g/L with the non validated cuvettes, which is not observed with STAGO cuvettes.

This variation shows that it is possible for a same sample to have normal result and pathologic result if non validated cuvettes are used, leading to a diagnosis risk for the patient.

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